

Anti-IL-6 Mouse Mab

Cat. No: DTA0003

DESCRIPTION

| | |
|-------------------------|--|
| Product Name | Anti-IL-6 Mouse Mab |
| Host Species | Mouse |
| Immunogen | Recombinant protein of human IL6 |
| Clonality | Monoclonal |
| Ig class | IgG |
| Purified by | Protein A/G affinity column |
| Applications | ELISA |
| Cross-Reactivity | Human |
| Conjugate | Unconjugated |
| Lot number | 9767804 |
| Concentration | 5.0mg/mL |
| Storage | For up to 36 months at -20°C. After receive, the antibody should be divided into the right amount. Avoid repeated freezing and thawing cycles. |
| Storage Buffer | PBS, pH 7.4 |

*This product is for research use only.