

Anti-IL-6 Rabbit Mab

Cat. No: DTA0004

DESCRIPTION

Product Name	Anti-IL-6 Rabbit Mab
Host Species	Rabbit
Immunogen	Recombinant protein of human IL6
Clonality	Monoclonal
Ig class	IgG
Purified by	Protein A/G affinity column
Applications	ELISA
Cross-Reactivity	Human
Conjugate	Unconjugated
Lot number	8628680
Concentration	5.0mg/mL
Storage	For up to 36 months at -20°C. After receive, the antibody should be divided into the right amount. Avoid repeated freezing and thawing cycles.
Storage Buffer	PBS, pH 7.4

*This product is for research use only.