

Anti-IFN- γ Mouse Mab

Cat. No: DTA0005

DESCRIPTION

| | |
|---------------------|---|
| Product Name | Anti-IFN-γ Mouse Mab |
|---------------------|---|

| | |
|---------------------|-------|
| Host Species | Mouse |
|---------------------|-------|

| | |
|------------------|------------|
| Clonality | Monoclonal |
|------------------|------------|

| | |
|-----------------|-----|
| Ig class | IgG |
|-----------------|-----|

| | |
|--------------------|-----------------------------|
| Purified by | Protein A/G affinity column |
|--------------------|-----------------------------|

| | |
|---------------------|-------|
| Applications | ELISA |
|---------------------|-------|

| | |
|-------------------------|-------|
| Cross-Reactivity | Human |
|-------------------------|-------|

| | |
|------------------|--------------|
| Conjugate | Unconjugated |
|------------------|--------------|

| | |
|-------------------|---------|
| Lot number | 5336361 |
|-------------------|---------|

| | |
|----------------------|----------|
| Concentration | 3.0mg/mL |
|----------------------|----------|

| | |
|----------------|--|
| Storage | For up to 36 months at -20°C. After receive, the antibody should be divided into the right amount. Avoid repeated freezing and thawing cycles. |
|----------------|--|

| | |
|-----------------------|-------------|
| Storage Buffer | PBS, pH 7.4 |
|-----------------------|-------------|

*This product is for research use only.