

# Anti-HLA-B27 Mouse Mab

Cat. No: DTA0018

## DESCRIPTION

<b>Product Name</b>	Anti-HLA-B27 Mouse Mab
<b>Host Species</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Ig class</b>	IgG
<b>Purified by</b>	Protein A/G affinity column
<b>Applications</b>	ELISA, FC
<b>Cross-Reactivity</b>	Human
<b>Conjugate</b>	Unconjugated/FITC
<b>Lot number</b>	7842474
<b>Concentration</b>	5.0mg/mL
<b>Storage</b>	For up to 36 months at -20°C. After receive, the antibody should be divided into the right amount. Avoid repeated freezing and thawing cycles.
<b>Storage Buffer</b>	PBS, pH 7.4

\*This product is for research use only.