

Anti-CD45 Mouse Mab

Cat. No: DTA0026

DESCRIPTION

Product Name	Anti-CD45 Mouse Mab
Host Species	Mouse
Clonality	Monoclonal
Ig class	IgG
Purified by	Protein A/G affinity column
Applications	ELISA, FC
Cross-Reactivity	Human
Conjugate	Unconjugated/FITC
Lot number	3860472
Concentration	3.0mg/mL
Storage	For up to 36 months at -20°C. After receive, the antibody should be divided into the right amount. Avoid repeated freezing and thawing cycles.
Storage Buffer	PBS, pH 7.4

*This product is for research use only.