

Guinea Pig IgG Control

Cat.No: DTA0167

DESCRIPTION

Name	Guinea Pig IgG Control
Purification	Protein A affinity column
Clonality	Polyclonal
Isotype	Guinea Pig IgG
Form	Liquid
Concentration	1 mg/ml
Purify	>90%
Storage Buffer	1×PBS, pH7.4, containing 1% BSA, 50% glycerol and 0.03% Proclin 300.
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.

Applications

Application	Recommended Usage
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!