

# Rabbit IgG[HRP] Control

Cat.No: DTA0182

## DESCRIPTION

Name	Rabbit IgG[HRP] Control
Conjugate	Horseradish peroxidase (HRP)
Clonality	Polyclonal
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 1% BSA,50% glycerol and 0.03% Proclin 300.
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.

## Applications

Application	Recommended Usage
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

**For Research Use Only!**