

# Rabbit Anti Protein C Tag Antibody[HRP], pAb

Cat.No: DTA0189

## DESCRIPTION

Name	Rabbit Anti Protein C Tag Antibody[HRP], pAb
Immunogen	Protein C tag peptide(EDQVDPRLIDGK) conjugated to keyhole limpet haemocyanin(KLH)
Conjugate	Horseradish peroxidase (HRP)
Host Species	Rabbit
Specificity	This antibody recognizes C-terminal and N-terminal region of Protein C tag fusion proteins.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 1%, 50% glycerol and 0.03% Proclin 300.
Clonality	Polyclonal
Isotype	Rabbit IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

## Applications

Application	Recommended Usage
ELISA	0.05-0.5μg/ml
Western Blot	0.2-1μg/ml
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

**For Research Use Only!**