

Rabbit Anti CBP Tag Antibody, pAb

Cat.No: DTA0208

DESCRIPTION

Name	Rabbit Anti CBP Tag Antibody, pAb
Immunogen	CBP Tag peptide(KRRWKKNFIAVSAANRFKKISSSGAL) conjugated to keyhole limpet haemocyanin(KLH)
Host Species	Rabbit
Specificity	This antibody recognizes CBP tag fusion proteins.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 50% glycerol and 0.03% Proclin 300.
Purification	Antigen affinity purification
Clonality	Polyclonal
Isotype	Rabbit IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	0.05-0.5 μg/ml
Western Blot	0.2-1 μg/ml
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!