

Mouse Anti DYKDDDDK Tag Antibody[HRP], mAb

Cat.No: DTA0217

DESCRIPTION

Name	Mouse Anti DYKDDDDK Tag Antibody[HRP], mAb
Immunogen	DYKDDDDK peptide conjugated to keyhole limpet haemocyanin(KLH)
Conjugate	Horseradish peroxidase (HRP)
Host Species	Mouse
Specificity	This antibody recognizes C-terminal, N-terminal, and internal region of DYKDDDDK tag fusion proteins.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 1% BSA, 50% glycerol and 0.03% Proclin 300.
Clonality	Monoclonal
Isotype	Mouse IgG2b, κ
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	0.1-1 μ g/ml
Western Blot	0.5-1 μ g/ml
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!