

# Mouse Anti GST Tag Antibody[Biotin], mAb

Cat.No: DTA0221

## DESCRIPTION

|                      |                                                                                                                                                                |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name                 | Mouse Anti GST Tag Antibody[Biotin], mAb                                                                                                                       |
| Immunogen            | Recombinant GST protein                                                                                                                                        |
| Conjugate            | Biotin                                                                                                                                                         |
| Host Species         | Mouse                                                                                                                                                          |
| Specificity          | This antibody recognizes C-terminal, N-terminal, and internal region of GST tag fusion proteins.                                                               |
| Form                 | Liquid                                                                                                                                                         |
| Concentration        | 0.5 mg/ml                                                                                                                                                      |
| Storage Buffer       | 1×PBS, pH7.4, containing 1% BSA, 50% glycerol and 0.03% Proclin 300.                                                                                           |
| Clonality            | Monoclonal                                                                                                                                                     |
| Isotype              | Mouse IgG1, $\kappa$                                                                                                                                           |
| Shipping and Storage | Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles. |
| Expiration Date      | Specified on product label.                                                                                                                                    |

## Applications

| Application        | Recommended Usage   |
|--------------------|---------------------|
| ELISA              | 0.05-0.5 $\mu$ g/ml |
| Western Blot       | 0.5-1 $\mu$ g/ml    |
| Other applications | User-optimized      |

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!