

# Mouse Anti V5 Tag Antibody[HRP], mAb

Cat.No: DTA0224

## DESCRIPTION

Name	Mouse Anti V5 Tag Antibody[HRP], mAb
Immunogen	V5 tag peptide conjugated to keyhole limpet haemocyanin(KLH)
Conjugate	Horseradish peroxidase (HRP)
Host Species	Mouse
Specificity	This antibody recognizes C-terminal and N-terminal region of V5 tag fusion proteins.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 0.1% BSA, 50% glycerol and 0.03% Proclin 300.
Purification	Purified by Protein A
Clonality	Monoclonal
Isotype	Mouse IgG1, $\kappa$
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

## Applications

Application	Recommended Usage
ELISA	0.02-0.2 $\mu$ g/ml
Western Blot	0.2-1 $\mu$ g/ml
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!