

# Rabbit Anti GFP Tag Antibody[HRP], pAb

Cat.No: DTA0227

## DESCRIPTION

Name	Rabbit Anti GFP Tag Antibody[HRP], pAb
Immunogen	recombinant GFP protein
Conjugate	Horseradish peroxidase (HRP)
Host Species	Rabbit
Specificity	This antibody recognizes GFP fusion proteins and reacts with various forms of GFP including EGFP, cGFP.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 1% BSA,50% glycerol and 0.03% Proclin 300.
Clonality	Polyclonal
Isotype	Rabbit IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

## Applications

Application	Recommended Usage
ELISA	0.05-0.5μg/ml
Western Blot	0.2-1μg/ml
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

**For Research Use Only!**