

Mouse Anti MMAE Antibody[HRP], mAb

Cat.No: DTA0241

DESCRIPTION

Name	Mouse Anti MMAE Antibody[HRP], mAb
Immunogen	MMAE conjugated to keyhole limpet haemocyanin(KLH)
Conjugate	Horseradish peroxidase (HRP)
Host Species	Mouse
Specificity	This antibody recognizes MMAE molecules and does not react with MMAF.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 0.1% BSA, 50% glycerol and 0.03% Proclin 300.
Clonality	Monoclonal
Isotype	Mouse IgG2a, κ
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	0.1-2 μ g/ml
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!