

Mouse Anti DXd Antibody[HRP], mAb

Cat.No: DTA0243

DESCRIPTION

Name	Mouse Anti DXd Antibody[HRP], mAb
Immunogen	DXd conjugated to keyhole limpet haemocyanin(KLH)
Conjugate	Horseradish peroxidase (HRP)
Host Species	Mouse
Specificity	This antibody recognizes DXd and DXd structurally similar molecules. The antibody can bind to DS-8201.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 0.05% BSA, 50% glycerol and 0.03% Proclin 300.
Purification	Purified by Protein A
Clonality	Monoclonal
Isotype	Mouse IgG1, λ
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	0.1-2 μ g/ml
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!