

# Mouse Anti MMAF Antibody[HRP], mAb

Cat.No: DTA0247

## DESCRIPTION

|                      |  |
|----------------------|--|
| Name                 | Mouse Anti MMAF Antibody[HRP], mAb   |
| Immunogen            | MMAF conjugated to keyhole limpet haemocyanin(KLH)   |
| Conjugate            | Horseradish peroxidase (HRP)   |
| Host Species         | Mouse  |
| Specificity          | This antibody recognizes MMAF molecules and does not react with MMAE.  |
| Form                 | Liquid   |
| Concentration        | 0.5 mg/ml  |
| Storage Buffer       | 1×PBS, pH7.4, containing 0.1% BSA, 50% glycerol and 0.03% Proclin 300.   |
| Clonality            | Monoclonal   |
| Isotype              | Mouse IgG1, $\lambda$  |
| Shipping and Storage | Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles. |
| Expiration Date      | Specified on product label.  |

## Applications

| Application        | Recommended Usage |
|--------------------|-------------------|
| ELISA              | 0.1-1 $\mu$ g/ml  |
| Other applications | User-optimized    |

感谢您选择德泰生物，我们竭诚为您服务！

**For Research Use Only!**