

Goat Anti-Rabbit IgG(H&L)[iFluor 488], pAb

Cat.No:DTA0255

DESCRIPTION

Name	Goat Anti-Rabbit IgG(H&L)[iFluor 488], pAb
Immunogen	Rabbit IgG (H&L)
Host Species	Goat
Specificity	This antibody can recognize Rabbit IgG.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 1% BSA, 0.03% Proclin 300 and 50% glycerol.
Conjugate	iFluor 488
Clonality	Polyclonal
Isotype	Goat IgG
Shipping and Storage	Shipped in ice bag. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
Flow Cytometry/ICC	1:200-1,500
IHC	1:500-1:2,500
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!