

Rabbit Anti-Rat IgG(H&L)[HRP], pAb

Cat.No:DTA0262

DESCRIPTION

Name	Rabbit Anti-Rat IgG(H&L)[HRP], pAb
Immunogen	Rat IgG (H&L)
Host Species	Rabbit
Specificity	This antibody can recognize Rat IgG.
Form	Liquid
Concentration	1 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 1% BSA, 0.03% Proclin 300 and 50% glycerol.
Conjugate	Horseradish peroxidase (HRP)
Clonality	Polyclonal
Isotype	Rabbit IgG
Shipping and Storage	Shipped in ice bag. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	1:10,000-1:50,000
Western Blot	1:5,000-1:20,000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!