

# Rabbit Anti-Rat IgG(H&L)[FITC], pAb

Cat.No:DTA0264

## DESCRIPTION

|                      |                                                                                                           |
|----------------------|-----------------------------------------------------------------------------------------------------------|
| Name                 | Rabbit Anti-Rat IgG(H&L)[FITC], pAb                                                                       |
| Immunogen            | Rat IgG (H&L)                                                                                             |
| Host Species         | Rabbit                                                                                                    |
| Specificity          | This antibody can recognize Rat IgG.                                                                      |
| Form                 | Liquid                                                                                                    |
| Concentration        | 1 mg/ml                                                                                                   |
| Storage Buffer       | 1×PBS, pH7.4, containing 1% BSA, 0.03% Proclin 300 and 50% glycerol.                                      |
| Conjugate            | FITC                                                                                                      |
| Clonality            | Polyclonal                                                                                                |
| Isotype              | Rabbit IgG                                                                                                |
| Shipping and Storage | Shipped in ice bag. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles. |
| Expiration Date      | Specified on product label.                                                                               |

## Applications

| Application        | Recommended Usage |
|--------------------|-------------------|
| Flow Cytometry/IF  | 1:200-1,500       |
| IHC                | 1:500-1:2,500     |
| Other applications | User-optimized    |

感谢您选择德泰生物，我们竭诚为您服务！

**For Research Use Only!**