

# Goat Anti-Human IgG(H&L)[Biotin], pAb

Cat.No:DTA0266

## DESCRIPTION

Name	Goat Anti-Human IgG(H&L)[Biotin], pAb
Immunogen	Human IgG (H&L)
Host Species	Goat
Specificity	This antibody can recognize Human IgG.
Form	Liquid
Concentration	1 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 1% BSA, 0.03% Proclin 300 and 50% glycerol.
Conjugate	Biotin
Clonality	Polyclonal
Isotype	Goat IgG
Shipping and Storage	Shipped in ice bag. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

## Applications

Application	Recommended Usage
ELISA	1:10,000-1:50,000
Western Blot	1:5,000-1:20,000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

**For Research Use Only!**