

Mouse Anti-Cat IgG Antibody, mAb

Cat.No: DTA0269

DESCRIPTION

Name	Mouse Anti-Cat IgG Antibody, mAb
Immunogen	Cat IgG
Host Species	Mouse
Specificity	This antibody reacts with the Cat IgG and has no cross-reactivity with mouse, rabbit, human, goat,dog,chicken and monkey immunoglobulins.
Form	Liquid
Concentration	0.5 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 0.1% BSA, 50% glycerol and 0.03% Proclin 300.
Clonality	Monoclonal
Isotype	Mouse IgG1, κ
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	1:2, 500–1:10, 000
WB	1:2, 500–1:10, 000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!