

Mouse Anti-Cat IgG Antibody[HRP], mAb

Cat.No: DTA0270

DESCRIPTION

Name	Mouse Anti-Cat IgG Antibody[HRP], mAb
Immunogen	Cat IgG (H&L)
Conjugate	Horseradish peroxidase (HRP)
Host Species	Mouse
Specificity	This antibody reacts with the Cat IgG and has no cross-reactivity with mouse, rabbit, human, goat,dog,chicken and monkey immunoglobulins.
Form	Liquid
Concentration	1 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 1% BSA, 0.03% Proclin 300 and 50% glycerol.
Clonality	Monoclonal
Isotype	Mouse IgG1, κ
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	1:20,000-1:50,000
Western Blot	1:20,000-1:40,000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!