

Mouse Anti-Chicken IgY Antibody[HRP], mAb

Cat.No: DTA0272

DESCRIPTION

Name	Mouse Anti-Chicken IgY Antibody[HRP], mAb
Immunogen	Chicken IgY(H&L)
Conjugate	Horseradish peroxidase (HRP)
Host Species	Mouse
Specificity	This antibody recognizes Chicken IgY.
Form	Liquid
Concentration	1 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 0.1% BSA, 50% glycerol and 0.03% Proclin 300.
Clonality	Monoclonal
Isotype	Mouse IgG1, κ
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	1:20,000–1:50,000
Western Blot	1:20,000–1:40,000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!