

Mouse Anti-Goat IgG Antibody, mAb

Cat.No: DTA0309

DESCRIPTION

Name	Mouse Anti-Goat IgG Antibody, mAb
Immunogen	Goat IgG (H&L)
Host Species	Mouse
Specificity	The antibody reacts with the Goat IgG and has no cross-reactivity with human, Rabbit, chicken, guinea pig, mouse and rat immunoglobulins.
Form	Liquid
Concentration	1 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 0.03% Proclin 300.
Purification	Purified by Protein A
Clonality	Monoclonal
Isotype	Mouse IgG1, κ
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8 °C. For long term storage, store it at -20 °C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	1:5,000–1:20,000
Western Blot	1:1,000–1:10,000
IP	1:1000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!