

Mouse Anti-Rabbit IgG Antibody, mAb

Cat.No: DTA0607

DESCRIPTION

Name	Mouse Anti-Rabbit IgG Antibody, mAb
Immunogen	Rabbit IgG (H&L)
Host Species	Mouse
Specificity	The antibody reacts with the rabbit IgG and has no cross-reactivity with human, goat, chicken, guinea pig, mouse, rat and monkey immunoglobulins. This antibody only reacts with native rabbit IgG and does not bind to the denatured or reduced rabbit IgG heavy chain or light chain.
Form	Liquid
Concentration	1 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 0.1% BSA, 50% glycerol and 0.03% Proclin 300.
Purification	Purified by Protein A
Clonality	Monoclonal
Isotype	Mouse IgG2a, κ
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2-8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	1:20,000-1:40,000
Western Blot	1:10,000-1:20,000
IP	1:1000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务！

For Research Use Only!