


## Anti CD19 Antibody, mAb

Cat. No: DTA0645

### Description

Product Name	Anti CD19 Antibody, mAb
Clone	SJ25C1
Specificity	This antibody recognizes CD19 molecules.
Form	Liquid
Storage Buffer	PBS, pH7.4
Purification	Purified by Protein A
Clonality	Monoclonal
Isotype	Mouse IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2 -8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

### Applications

Application	Recommended Usage
Flow Cyt	
Other applications	User -optimized

Note: We can offer this product according to your requirement, including product size, buffer components, etc.

*For research use only. Not intended for human and animal therapeutic or diagnostic use.*