

## Anti CD45RO Antibody, mAb

Cat. No: DTA0655

### Description

|                      |                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti CD45RO Antibody, mAb                                                                                                                                     |
| Clone                | UCHL1                                                                                                                                                         |
| Specificity          | This antibody recognizes CD45RO molecules.                                                                                                                    |
| Form                 | Liquid                                                                                                                                                        |
| Storage Buffer       | PBS, pH7.4                                                                                                                                                    |
| Purification         | Purified by Protein A                                                                                                                                         |
| Clonality            | Monoclonal                                                                                                                                                    |
| Isotype              | Mouse IgG                                                                                                                                                     |
| Shipping and Storage | Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2 -8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles. |
| Expiration Date      | Specified on product label.                                                                                                                                   |

### Applications

| Application        | Recommended Usage |
|--------------------|-------------------|
| Flow Cyt           | √                 |
| IHC                | √                 |
| Other applications | User -optimized   |

Note: We can offer this product according to your requirement, including product size, buffer components, etc.

---

*For research use only. Not intended for human and animal therapeutic or diagnostic use.*