

Anti hOX40L Antibody, mAb

Cat. No: DTA0664

Description

Product Name	Anti hOX40L Antibody, mAb
Clone	159403
Specificity	This antibody recognizes hOX40L molecules.
Form	Liquid
Storage Buffer	PBS, pH7.4
Purification	Purified by Protein A
Clonality	Monoclonal
Isotype	Mouse IgG
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2 -8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
Flow Cyt	√
Other applications	User -optimized

Note: We can offer this product according to your requirement, including product size, buffer components, etc.

For research use only. Not intended for human and animal therapeutic or diagnostic use.