

## Anti IL12 Antibody, mAb

Cat. No: DTA0673

### Description

Product Name	Anti IL12 Antibody, mAb
Recommended Matches	DTA0674
Specificity	This antibody recognizes IL12 molecules.
Form	Liquid
Storage Buffer	PBS, pH7.4
Purification	Purified by Protein A
Clonality	Monoclonal
Isotype	Mouse IgG1
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2 -8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

### Applications

Application	Recommend
Flow Cyt、ELISA、IHC	√
Other applications	User -optimized

Note: We can offer this product according to your requirement, including product size, buffer components, etc.

*For research use only. Not intended for human and animal therapeutic or diagnostic use.*