

## Anti IFN- $\gamma$ Antibody, mAb

Cat. No: DTA0675

### Description

|                      |                                                                                                                                                               |
|----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti IFN- $\gamma$ Antibody, mAb                                                                                                                              |
| Specificity          | This antibody recognizes IFN- $\gamma$ molecules.                                                                                                             |
| Form                 | Liquid                                                                                                                                                        |
| Storage Buffer       | PBS, pH7.4                                                                                                                                                    |
| Purification         | Purified by Protein A                                                                                                                                         |
| Clonality            | Monoclonal                                                                                                                                                    |
| Isotype              | Mouse IgG1                                                                                                                                                    |
| Shipping and Storage | Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2 -8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles. |
| Expiration Date      | Specified on product label.                                                                                                                                   |

### Applications

| Application        | Recommend       |
|--------------------|-----------------|
| Flow Cyt、ELISA     | √               |
| Other applications | User -optimized |

Note: We can offer this product according to your requirement, including product size, buffer components, etc.

*For research use only. Not intended for human and animal therapeutic or diagnostic use.*