

Anti TNF- α Antibody, mAb

Cat. No: DTA0678

Description

Product Name	Anti TNF- α Antibody, mAb
Recommended Matches	DTA0677
Specificity	This antibody recognizes TNF- α molecules.
Form	Liquid
Storage Buffer	PBS, pH7.4
Purification	Purified by Protein A
Clonality	Monoclonal
Isotype	Mouse IgG1
Shipping and Storage	Shipped in ice bag. The antibody can be stored for 1-2 weeks at 2 -8°C. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommend
Flow Cyt、ELISA、IHC	√
Other applications	User -optimized

Note: We can offer this product according to your requirement, including product size, buffer components, etc.

For research use only. Not intended for human and animal therapeutic or diagnostic use.