

Recombinant Human Serum Albumin (rHSA)

Cat.No:DTD0300

DESCRIPTION

Name Recombinant Human Serum Albumin (rHSA)

Molecular weight 66.4 kDa

Purity $\geq 99.0\%$

Expression System Pig

Endotoxins $< 0.01\text{EU/mg}$



Residual Exogenous DNA $\leq 0.5\text{ng/g}$

pH 6.8

Physical Appearance Off-white to white powder.

Stability & Storage Store at 2°C to 8°C. Valid for 2 years.

LC-MS



Thank you for choosing Detai Biological. We are dedicated to serving you!

For Research Use Only