Streptavidin-HRP

Cat.No: DTA0297

DESCRIPTION

Name	Streptavidin-HRP
Size	1 mg
Conjugate	Horseradish peroxidase (HRP)
Form	Liquid
Concentration	1 mg/ml
Storage Buffer	1×PBS, pH7.4, containing 1% BSA, 50% glycerol and 0.03% Proclin 300.
Shipping and Storage	Shipped in ice bag. Store in the dark. For long term storage, store it at -20°C. Avoid repeated freezing and thawing cycles.
Expiration Date	Specified on product label.

Applications

Application	Recommended Usage
ELISA	1:5,000-1:20,000
Western Blot/Dot Blot	1:2,000-1:20,000
Other applications	User-optimized

感谢您选择德泰生物，我们竭诚为您服务!

For Research Use Only!